

# Developing and staining slides for mRNA detection in mouse embryo tissue sections

## Antonio Simeone

MRC Centre for Developmental Neurobiology, New Hunt's House, King's College London, Guy's Campus, London Bridge, London SE1 1UL UK.

International Institute of Genetics and Biophysics, CNR, Via G. Marconi, 12, 80125 Naples, Italy.

### Reagents

- ◆ Developing solution: Kodak D19 developer at 160 g/litre
- ◆ Fixing solution: Kodak Unifix at 130 g/litre

### A. Developing slides

- 1 Remove box with slides from 4 °C and leave 1 h at room temperature.
- 2 In the dark-room transfer the slides to a slide rack.<sup>a</sup>
- 3 Place the slide rack in a freshly prepared developing solution for 2 min, in 1% glycerol:1% acetic acid (alternatively a wash in water is sufficient) for 1 min, then in fixing solution for 5–10 min, all at room temperature.
- 4 Wash slides twice in a large volume of water for 10 min each.
- 5 Place slides in 0.02% toluidine blue solution<sup>b</sup> for 5 min, then wash in water for 1–2 min.
- 6 Dehydrate and destain<sup>c</sup> slides by passing through 30%, 60%, 85%, 95% ethanol, then twice in 100% ethanol, and then air dry.
- 7 Mount the slides with a clean coverslip using a few drops of a mounting agent (Permount or DPX) and dry in a fume-hood.

### B. Rescuing slides from high background in the emulsion

**If a high level of background is detected in the emulsion upon checking one or two experimental slides, proceed as follows with the remaining slides.**

1. Transfer the undeveloped slides to a slide rack.
2. Place the slide rack in a low concentration fixing solution (30 g/litre) for 20 min; the slides will appear completely transparent.
3. Wash four times in a large volume of 0.1 × SSC for 20 min each.

4. Dehydrate slides through 30%, 60%, 85%, 95%, then twice in 100% ethanol, and then air dry.
5. Repeat autoradiography procedure ( see Autoradiography for mRNA detection in mouse embryo tissue sections<sup>a</sup>) and develop again.

### Notes

- a Before developing all the slides, check one or two of them to verify that there is no background.
- b Alternative stains can be used such as eosin, haematoxylin/eosin, or Giemsa.
- c To better visualize faint signals we suggest a pale blue stain.