

Subbing of slides for mouse embryo sections for mRNA detection

Antonio Simeone

MRC Centre for Developmental Neurobiology, New Hunt's House, King's College London, Guy's Campus, London Bridge, London SE1 1UL UK.
International Institute of Genetics and Biophysics, CNR, Via G. Marconi, 12, 80125 Naples, Italy.

Equipment and reagents

- ◆ Slide racks for 25 slides
- ◆ Oven
- ◆ TESPA (aminopropyltriethoxysilane) (Sigma, A3648)
- ◆ Chrome-alum gelatin: 0.5% chrome-alum, 0.5% gelatin. Dissolve gelatin at a concentration of 5 g per litre in DEPC treated water at 60 °C on a stirrer, and then add 5 g chrome-alum per litre. Filter through Whatman 3MM paper.

Method A

- 1 Put slides in slide rack(s).
- 2 Dip for 30 sec in 10% HCl, 70% ethanol.
- 3 Wash in distilled water.
- 4 Dip for 30–60 sec in 95% ethanol.
- 5 Dry in an oven at 70–80 °C for 20–40 min and allow to cool.
- 6 Dip slides in 2% TESPA in acetone for 30 sec.
- 7 Wash twice in acetone and then in distilled water.
- 8 Dry at 37–42 °C for several hours. Can store airtight for months at room temperature.

Method B

1. Wash slides with an abundant volume of 0.5% SDS.
2. Wash slides in tap-water for 1 h.
3. Put slides in 25 slide rack(s) and wash twice in bidistilled water.
4. Dip slides in chrome-alum gelatin for 2 min.
5. Air dry and bake at 150 °C for 2–3 h or at 37 °C overnight. Can store airtight for months at room temperature.